

Suprapubic & Retropubic Prostatectomy Set

14 PCS (RAS-113-13)



Fig. 01
RG-270-20



Fig. 02
RI-310-02



Fig. 03
RAD-122-00



Fig. 04
RI-350-04



Fig. 05
RE-122-23



Fig. 06
RE-123-23



Fig. 25
RAS-113-13
Stainless Steel Box

www.rheingroup.com

812

Rhein Group

Sr #	Qty	Description
Fig-01	06-Pcs	Rochester-Pean Forceps, Curved 20cm / 8"
Fig-02	02-Pcs	Richardson Retractor, Blade, 36 x 28 mm 24cm / 9 1/2"
Fig-03	01-Pcs	Masson-Judd Bladder Retractor, with Center Blade 25x70 mm
Fig-04	02-Pcs	Kelly Retractor, Blade, 50 x 35 mm 26cm / 10 1/4"
Fig-05	01-Pcs	Mayo Dissecting Scissors, Straight, TC Blades 23cm / 9"
Fig-06	01-Pcs	Mayo Dissecting Scissors, Curved, TC Blades 23cm / 9"
Fig-07	01-Pcs	Stainless Steel Instrument Box