

Laparoscopy Set

21 PCS (RAS-122-22)



Fig. 01
RD-122-03



Fig. 02
RF-250-12



Fig. 03
RC-204-15



Fig. 04
RC-586-13



Fig. 05
RS-288-01



Fig. 06
RAE-116-02



Fig. 07
RAE-236-01



Fig. 08
RAE-336-32



Fig. 09
RAE-390-25



Fig. 10
RH-154-25



Fig. 11
RAA-104-15



Fig. 12
RE-122-17



Fig. 13
RL-130-18



Fig. 14
RAS-122-22
Stainless Steel Box

Laparoscopy Set

21 PCS (RAS-122-22)

Sr #	Qty	Description
Fig-01	01-Pcs	Standard Scalpal Handle, Solid, Straight No.3
Fig-02	01-Pcs	Adson Tissue Forceps, 1 x 2 Teeth 12cm / 4 3/4"
Fig-03	02-Pcs	Leriche Artery Forceps, Curved 15cm / 6"
Fig-04	06-Pcs	Backhaus Towel Clamp Forceps 13cm / 5 1/4"
Fig-05	01-Pcs	Heaney Retractor, 82 x 28 mm 27cm / 10 3/4"
Fig-06	01-Pcs	Graves Vaginal Speculum, Medium Size, 95 x 35 mm
Fig-07	01-Pcs	Auvard Weighted Speculum, Light 21lbs
Fig-08	01-Pcs	Sims Uterine Sound Malleable 32cm / 12 1/2"
Fig-09	01-Pcs	Schroeder Tenaculum Forceps, Straight 25cm / 10"
Fig-10	01-Pcs	Foerster Sponge Forceps, Straight, Serrated 25cm / 10"
Fig-11	02-Pcs	Allis Tissue Forceps, 4 x 5 Teeth 15cm / 6"
Fig-12	01-Pcs	Mayo Dissecting Scissors, Straight, TC Blades 17cm / 6 3/4"
Fig-13	01-Pcs	Crile-Wood, Needle Holder, TC Jaws 18cm / 7"
Fig-14	01-Pcs	Stainless Steel Instrument Box