

# Basic Craniotomy Set

104 PCS (RAS-125-25)



Fig. 31  
RP-156-19



Fig. 32  
RP-152-19



Fig. 33  
RP-248-17



Fig. 34  
RP-186-20



Fig. 35  
RP-230-15



Fig. 36  
RP-172-13



Fig. 37  
RP-318-03



Fig. 38  
RP-330-03



Fig. 39  
RP-336-03



Fig. 40  
RP-194-20



Fig. 41  
RW-316-20



Fig. 42  
RAA-104-15



Fig. 43  
RP-114-50



Fig. 44  
RP-112-00



Fig. 45  
RN-128-27

# Basic Craniotomy Set

104 PCS (RAS-125-25)



Fig. 46  
RN-128-00



Fig. 47  
RN-130-02



Fig. 48  
RO-396-01



Fig. 49  
RO-330-04



Fig. 50  
RO-338-01



Fig. 51  
RO-634-17



Fig. 52  
RO-658-18



Fig. 53  
RO-690-23



Fig. 54  
RE-122-14



Fig. 55  
RE-123-14



Fig. 56  
RE-203-18



Fig. 57  
RL-130-15



Fig. 58  
RL-130-18



Fig. 59  
RAS-125-25  
Stainless Steel Box

# Basic Craniotomy Set

104 PCS (RAS-125-25)

| Sr #   | Qty    | Description                                                    |
|--------|--------|----------------------------------------------------------------|
| Fig-01 | 01-Pcs | Frazier Suction Tube, 2.3 mm Fr.7 19cm / 7 1/2"                |
| Fig-02 | 01-Pcs | Frazier Suction Tube, 3 mm / Fr.9 19cm / 7 1/2"                |
| Fig-03 | 01-Pcs | Scalpel Handle, Solid, Straight No.3                           |
| Fig-04 | 01-Pcs | Adson Dressing Forceps, Serrated 12cm / 4 3/4"                 |
| Fig-05 | 01-Pcs | Adson Tissue Forceps, 1 x 2 Teeth 12cm / 4 3/4"                |
| Fig-06 | 01-Pcs | Semken Tissue Forceps, 1 x 2 Teeth 15cm / 6"                   |
| Fig-07 | 01-Pcs | Cushing Tissue Forceps, Straight, 1 x 2 Teeth 17cm / 6 3/4"    |
| Fig-08 | 01-Pcs | Taylor Dressing Forceps, Serrated 18.5cm / 7 1/4"              |
| Fig-09 | 06-Pcs | Halsted Mosquito Forceps, Straight, Delicate 12cm / 4 3/4"     |
| Fig-10 | 06-Pcs | Halsted Mosquito Forceps, Curved, Delicate 12cm / 4 3/4"       |
| Fig-11 | 12-Pcs | Dandy Artery Forceps, Curved Side ways 14cm / 5 1/2"           |
| Fig-12 | 02-Pcs | Rochester-Ochsner Forceps, Straight, 1 x 2 Teeth 16cm / 6 1/4" |
| Fig-13 | 02-Pcs | Rochester-Ochsner Forceps, Straight, 1 x 2 Teeth 18cm / 7"     |
| Fig-14 | 02-Pcs | Adson Artery Forceps, Straight 18cm / 7"                       |
| Fig-15 | 02-Pcs | Adson Artery Forceps, Curved 18cm / 7"                         |
| Fig-16 | 08-Pcs | Backhaus Towel Clamp Forceps 9cm / 3 1/2"                      |
| Fig-17 | 08-Pcs | Backhaus Towel Clamp Forceps 13cm / 5"                         |
| Fig-18 | 02-Pcs | Senn Miller Retractor, Sharp, Double ended 16cm / 6 1/4"       |
| Fig-19 | 02-Pcs | Muphry Retractor, 4 Prong, Sharp 19cm / 7 1/2"                 |
| Fig-20 | 01-Pcs | Jansen-Wagner Retractor, Sharp 10cm / 4"                       |
| Fig-21 | 01-Pcs | Jansen Retractor, Blunt 10cm / 4"                              |
| Fig-22 | 01-Pcs | Adson-Baby Retractor, Sharp, 3 x 4 Teeth 20cm / 8"             |
| Fig-23 | 01-Pcs | Beckman Retractor, Sharp 31cm / 12 1/4"                        |
| Fig-24 | 01-Pcs | Foerster Sponge Forceps, Straight, Serrated 25cm / 10"         |
| Fig-25 | 01-Pcs | Kelly Scissors, Straight, Sharp Points 16cm / 6 1/4"           |
| Fig-26 | 01-Pcs | Potts-Smith Scissors, 25 degree Angled 19cm / 7 1/2"           |
| Fig-27 | 01-Pcs | Universal Wire Cutting Scissors, Serrated 12cm / 4 3/4"        |
| Fig-28 | 01-Pcs | Martel Saw Guide 33cm / 13"                                    |
| Fig-29 | 01-Pcs | Love Nerve Root Retractor, 45 degree Angled 19cm / 7 3/4"      |
| Fig-30 | 01-Pcs | Scoville Brain Spatula 20cm / 8"                               |
| Fig-31 | 01-Pcs | Love Nerve Root Retractor, 90 degree Angled 19cm / 7 3/4"      |
| Fig-32 | 01-Pcs | Love Nerve Root Retractor, Straight 19cm / 7 3/4"              |
| Fig-33 | 01-Pcs | Cushing Spatula Spoon 17cm / 6 3/4"                            |
| Fig-34 | 01-Pcs | Dandy Nerve Hook, Straight 20cm / 8"                           |
| Fig-35 | 01-Pcs | Frazier Dura Dissector 15cm / 6"                               |
| Fig-36 | 02-Pcs | Frazier Dura Hook, Sharp 13cm / 5"                             |
| Fig-37 | 01-Pcs | Colclough Rongeur, 3 mm Upwards Bite 15cm / 6"                 |
| Fig-38 | 01-Pcs | Ferris-Smith-Kerrison Cervical Rongeur, 3 mm Cutting 18cm / 7" |
| Fig-39 | 01-Pcs | Cushing Laminectomy Rongeur, 2 x 10 mm, Straight 18cm / 7"     |
| Fig-40 | 01-Pcs | Adson Dura Hook, Sharp 20cm / 8"                               |
| Fig-41 | 01-Pcs | Freer Elevator, Double End, Sharp / Blunt 16cm / 6 1/4"        |
| Fig-42 | 02-Pcs | Allis Tissue Forceps, 4 x 5 Teeth 15cm / 6"                    |
| Fig-43 | 02-Pcs | Gigli Wire Saw 50cm / 20"                                      |
| Fig-44 | 02-Pcs | Gigli Saw Handle                                               |
| Fig-45 | 01-Pcs | Hudson Brace 27cm / 10 3/4"                                    |
| Fig-46 | 01-Pcs | Hudson Extension Piece                                         |
| Fig-47 | 01-Pcs | Drill Burr 14mm                                                |
| Fig-48 | 01-Pcs | Adson Periosteal Elevator, Straight, Sharp 17cm / 6 3/4"       |
| Fig-49 | 01-Pcs | Cushing Periosteal Elevator, 4 mm 17cm / 6 3/4"                |
| Fig-50 | 01-Pcs | Langenbeck Periosteal Elevator, 16 mm 19cm / 7 1/2"            |
| Fig-51 | 01-Pcs | Luer Bone Rongeur Forceps, Straight 17cm / 6 3/4"              |
| Fig-52 | 01-Pcs | Zaufal-Jansen Bone Rongeur Forceps 18cm / 7"                   |
| Fig-53 | 01-Pcs | Stille Bone Rongeur Forceps, Angled 23cm / 9"                  |
| Fig-54 | 01-Pcs | Mayo Dissecting Scissors, Straight, TC Blades 14.5cm / 5 3/4"  |
| Fig-55 | 01-Pcs | Mayo Dissecting Scissors, Curved, TC Blades 14.5cm / 5 3/4"    |
| Fig-56 | 01-Pcs | Metzenbaum Scissors, Curved, TC Blades 18cm / 7"               |
| Fig-57 | 01-Pcs | Crile-Wood, Needle Holder, TC Jaws 15cm / 6"                   |
| Fig-58 | 01-Pcs | Crile-Wood, Needle Holder, TC Jaws 18cm / 7"                   |
| Fig-59 | 01-Pcs | Stainless Steel Instrument Box                                 |