


Basic Laparotomy Set

110 PCS (RAS-101-01)




Basic Laparotomy Set

110 PCS (RAS-101-01)



Basic Laparotomy Set

110 PCS (RAS-101-01)



Basic Laparotomy Set

110 PCS (RAS-101-01)

Sr #	Qty	Description
Fig-01	01-Pcs	Yankauer Suction Tube 27cm / 11"
Fig-02	01-Pcs	Poole Suction Tube, Straight, Diameter 10mm / 3/32"
Fig-03	02-Pcs	Scalpel Handle, Solid, Straight No.3
Fig-04	02-Pcs	Scalpel Handle, Solid, Straight No.4
Fig-05	01-Pcs	Scalpel Handle, Solid, Straight No.7
Fig-06	01-Pcs	Dressing Forceps, Standard 14.5cm 5 1/2"
Fig-07	01-Pcs	Tissue Forceps, 1x2 Teeth 14.5cm 5 1/2"
Fig-08	01-Pcs	Adson Dressing Forceps, Serrated 12cm / 4 3/4"
Fig-09	01-Pcs	Adson Tissue Forceps, 1x2 Teeth 12cm / 4 3/4"
Fig-10	01-Pcs	Cushing Dressing Forceps, Serrated, Straight 17.5cm / 7"
Fig-11	01-Pcs	Cushing Dressing Forceps, 1x2 Teeth, Straight 17cm / 6 3/4"
Fig-12	01-Pcs	Russ-Model Dissecting Forceps 20cm / 8"
Fig-13	06-Pcs	Halsted Mosquito Forceps, Straight 12.5cm 5"
Fig-14	06-Pcs	Halsted Mosquito Forceps, Curved 12.5cm 5"
Fig-15	06-Pcs	Kelly Forceps, Straight 14cm / 5 1/2"
Fig-16	12-Pcs	Kelly Forceps, Curved 14cm / 5 1/2"
Fig-17	06-Pcs	Rochester-Pean Forceps, Curved 16cm / 6 1/4"
Fig-18	02-Pcs	Rochester-Pean Forceps, Curved 20cm / 8"
Fig-19	04-Pcs	Rochester-Ochsner Forceps, Straight, 1x2 Teeth 16cm / 6 1/4"
Fig-20	01-Pcs	Baby Mixer Forceps, Full Curved 18cm / 7"
Fig-21	01-Pcs	Mixer Forceps 23cm / 9"
Fig-22	01-Pcs	Mixer Forceps, Full Curved 18cm / 7"
Fig-23	02-Pcs	Lahey Gall Duct Forceps 19cm / 7 1/2"
Fig-24	04-Pcs	Schmidt Tonsil Forceps, Slight Curved 19cm / 7 1/2"
Fig-25	08-Pcs	Backhaus Towel Clamp Forceps 13cm / 5 1/4"
Fig-26	01-Pcs	Parker-Langenbeck Retractor, Set of 2 21cm / 8 1/4"
Fig-27	02-Pcs	Goel Retractor, 40x 40 mm, 30 x 30 mm 19cm / 7 1/2"
Fig-28	01-Pcs	Deaver Retractor, 25mm 30cm / 12"
Fig-29	01-Pcs	Deaver Retractor, 50mm 30cm / 12"
Fig-30	01-Pcs	Richardson Retractor, Blade, 28 x 20 mm 24cm / 9 1/2"
Fig-31	01-Pcs	Richardson Retractor, Blade, 44 x 38 mm 24cm / 9 1/2"
Fig-32	01-Pcs	Kelly Retractor, Blade, 65 x 50 mm 26cm / 10 1/4"
Fig-33	01-Pcs	Balfour Retractor With Interchangeable, Blades 180 mm
Fig-34	01-Pcs	Foenster Sponge Forceps, Straight, Serrated 25cm / 10"
Fig-35	01-Pcs	Operating Scissors, Straight, Sharp/Blunt 14.5cm 5 1/2"
Fig-36	01-Pcs	Ribbon Retractor, Flexible, 30 mm 33cm / 13"
Fig-37	01-Pcs	Ribbon Retractor, Flexible, 40 mm 33cm / 13"
Fig-38	01-Pcs	Ribbon Retractor, Flexible, 50 mm 33cm / 13"
Fig-39	01-Pcs	Forceps Holder
Fig-40	01-Pcs	Forceps Holder
Fig-41	01-Pcs	Forceps Holder
Fig-42	01-Pcs	DeBakey Tissue Forceps 20cm / 8"
Fig-43	02-Pcs	Babcock Tissue Forceps, Straight 16cm / 6 1/4"
Fig-44	02-Pcs	Babcock Tissue Forceps, Straight 24cm / 9 1/2"
Fig-45	02-Pcs	Allis Tissue Forceps, 4 x 5 Teeth 15cm / 6"
Fig-46	02-Pcs	Allis Tissue Forceps, 5 x 6 Teeth 25cm / 10"
Fig-47	01-Pcs	Mayo Operating Scissors, Straight, TC Blades 17cm / 6 3/4"
Fig-48	01-Pcs	Mayo Operating Scissors, Curved, TC Blades 17cm / 6 3/4"
Fig-49	01-Pcs	Metzenbaum Scissors, Straight, TC Blades 23cm / 9"
Fig-50	01-Pcs	Metzenbaum Scissors, Curved, TC Blades 18cm / 7"
Fig-51	01-Pcs	Metzenbaum Scissors, Curved, TC Blades 23cm / 9"
Fig-52	01-Pcs	Mayo-Hegar, Needle Holder, TC Jaws 16cm / 6 1/4"
Fig-53	01-Pcs	Mayo-Hegar, Needle Holder, TC Jaws 18cm / 7"
Fig-54	01-Pcs	Mayo-Hegar, Needle Holder, TC Jaws 20cm / 8"
Fig-55	01-Pcs	Mayo-Hegar, Needle Holder, TC Jaws 26cm / 10 1/4"
Fig-56	01-Pcs	Stainless Steel Ruler 15cm / 6"
Fig-57	01-Pcs	Stainless Steel Instrument Box