

Amputation Set

16 PCS (RAS-133-33)



Fig. 01
RD-154-20



Fig. 02
RP-110-33



Fig. 03
RO-306-27



Fig. 04
RP-114-50



Fig. 05
RP-112-00



Fig. 06
RN-104-31



Fig. 07
RL-144-06



Fig. 08
RL-144-13



Fig. 09
RL-144-19



Fig. 10
RL-144-25



Fig. 11
RL-144-32



Fig. 12
RO-338-01



Fig. 13
RO-634-17



Fig. 14
RO-722-19



Fig. 15
RAS-133-33
Stainless Steel Box

Amputation Set

16 PCS (RAS-133-33)

Sr #	Qty	Description
Fig-01	01-Pcs	Liston Amputation Knife, 20cm / 8"
Fig-02	01-Pcs	Martel Gigli Saw Guide, 33cm / 13"
Fig-03	01-Pcs	Putti Bone Rasp, Flat Ends, 27cm / 10 1/2"
Fig-04	02-Pcs	Gigli Wire Saw, 50cm / 20"
Fig-05	01-Pcs	Gigli Strully Saw Handle
Fig-06	01-Pcs	Satterlee Saw, Blade 23cm / 9" 30cm / 12"
Fig-07	01-Pcs	Hibbs / Hilder Osteotome, Straight, 6mm 24cm / 9 1/2"
Fig-08	01-Pcs	Hibbs / Hilder Osteotome, Straight, 13mm 24cm / 9 1/2"
Fig-09	01-Pcs	Hibbs / Hilder Osteotome, Straight, 19mm 24cm / 9 1/2"
Fig-10	01-Pcs	Hibbs / Hilder Osteotome, Straight, 25mm 24cm / 9 1/2"
Fig-11	01-Pcs	Hibbs / Hilder Osteotome, Straight, 32mm 24cm / 9 1/2"
Fig-12	01-Pcs	Langenbeck Periosteal Elevator, Wide Sharp 16mm 19cm / 7 1/2"
Fig-13	01-Pcs	Luer Bone Rongeur Forceps, Straight 17cm / 6 3/4"
Fig-14	01-Pcs	Liston Bone Cutting Forceps, Straight 19cm / 7 1/2"
Fig-15	01-Pcs	Stainless Steel Instrument Box