

Closed Thoracostomy Set

14 PCS (RAS-136-36)



Fig. 01
RC-114-25



Fig. 02
RC-114-30



Fig. 03
RD-122-03




Fig. 04
RC-176-16



Fig. 05
RC-586-13



Fig. 06
RH-134-25




Fig. 07
RE-122-14



Fig. 08
RAS-136-36
Stainless Steel Box

Sr #	Qty	Description
Fig-01	01-Pcs	Nelson Trocar 21cm / 8 1/4"
Fig-02	01-Pcs	Nelson Trocar 21cm / 8 1/4"
Fig-03	01-Pcs	Standard Scalpel Handle, Solid, Straight, No. 3
Fig-04	02-Pcs	Crile-Rankin Forceps, Straight 16cm / 6 1/4"
Fig-05	06-Pcs	Backhaus Towel Clamp Forceps 13cm / 5 1/4"
Fig-06	01-Pcs	Foerster Sponge Forceps, Straight, Serrated 25cm / 10"
Fig-07	01-Pcs	Mayo Dissecting Scissors, Straight, TC Blades 14.5cm / 5 3/4"
Fig-08	01-Pcs	Stainless Steel Instrument Box