

# Basic Open Heart Set

93 PCS (RAS-139-39)



Fig. 01  
RC-250-27

Fig. 02  
RD-122-03

Fig. 03  
RD-132-04

Fig. 04  
RD-142-07

Fig. 05  
RF-248-12



Fig. 06  
RF-234-17

Fig. 07  
RC-166-14

Fig. 08  
RC-204-15

Fig. 09  
RC-232-16

Fig. 10  
RC-138-12



Fig. 11  
RC-270-22

Fig. 12  
RC-272-20

Fig. 13  
RC-346-23

Fig. 14  
RC-324-18

Fig. 15  
RC-586-13

# Basic Open Heart Set

93 PCS (RAS-139-39)



Fig. 16  
RC-564-20

Fig. 17  
RC-336-20

Fig. 18  
RS-314-00 Set

Fig. 19  
RH-154-25

Fig. 20  
RE-195-14



Fig. 21  
RE-211-28

Fig. 22  
RE-211-30

Fig. 23  
RR-116-20

Fig. 24  
RR-118-16

Fig. 25  
RR-276-22




Fig. 26  
RR-276-21

Fig. 27  
RR-324-20


Fig. 28  
RR-324-25

Fig. 29  
RR-274-25

Fig. 30  
RAA-106-25

# Basic Open Heart Set

93 PCS (RAS-139-39)



# Basic Open Heart Set

93 PCS (RAS-139-39)

| Sr #   | Qty    | Description                                                         |
|--------|--------|---------------------------------------------------------------------|
| Fig-01 | 02-Pcs | Yankauer Suction Tube, Chrome Plated 23cm / 9"                      |
| Fig-02 | 01-Pcs | Standard Scalpal Handle, Solid, Straight No. 3                      |
| Fig-03 | 01-Pcs | Standard Scalpal Handle, Solid, Straight No. 4                      |
| Fig-04 | 01-Pcs | Scalpel Handle, Solid, Straight No. 7                               |
| Fig-05 | 02-Pcs | Adson Tissue Forceps, 1 x 2 Teeth 12cm / 4 3/4"                     |
| Fig-06 | 01-Pcs | Cushing Tissue Forceps, Straight, 1 x 2 Teeth 17cm / 6 3/4"         |
| Fig-07 | 06-Pcs | Crile Forceps, Curved 14cm / 5 1/2"                                 |
| Fig-08 | 04-Pcs | Lerich Artery Forceps, Curved 15cm / 6"                             |
| Fig-09 | 06-Pcs | Pean Artery Forceps, Straight 16cm / 6 1/4"                         |
| Fig-10 | 18-Pcs | Halsted Mosquito Forceps, Curved 12.5cm / 5"                        |
| Fig-11 | 04-Pcs | Rochester-Pean Forceps, Curved 22cm / 8 3/4"                        |
| Fig-12 | 02-Pcs | Rochester Ochsner Forceps, b2 Teeth, Straight 20cm / 8"             |
| Fig-13 | 01-Pcs | Gemini-Mixer Forceps 23cm / 9"                                      |
| Fig-14 | 02-Pcs | Baby-Mixer Forceps, Full Curved 18cm / 7"                           |
| Fig-15 | 08-Pcs | Backhaus Towel Clamp Forceps 13cm / 5 1/4"                          |
| Fig-16 | 04-Pcs | Tubing Clamps, With guard 20cm / 8"                                 |
| Fig-17 | 02-Pcs | Heiss Forceps, Curved 20cm / 8"                                     |
| Fig-18 | 01-Pcs | Richardson-Eastman Retractor, Set of 2, Double ended 26cm / 10 1/4" |
| Fig-19 | 02-Pcs | Foerster Sponge Forceps, Straight, Serrated 25cm / 10"              |
| Fig-20 | 01-Pcs | Metzenbaum Scissors, Curved, TC Blades 14cm / 5 1/2"                |
| Fig-21 | 01-Pcs | Metzenbaum Scissors, Curved, TC Blades 28cm / 11"                   |
| Fig-22 | 01-Pcs | Metzenbaum Scissors, Curved, TC Blades 30cm / 12"                   |
| Fig-23 | 01-Pcs | DeBakey Thoracic Tissue Forceps 20cm / 8"                           |
| Fig-24 | 01-Pcs | DeBakey Tissue Forceps 16cm / 6 1/4"                                |
| Fig-25 | 01-Pcs | Potts Cardiovascular Forceps 22cm / 8 3/4"                          |
| Fig-26 | 02-Pcs | Potts Cardiovascular Forceps, Straight 21.5cm / 8 1/2"              |
| Fig-27 | 01-Pcs | DeBakey Cardiovascular Forceps 20cm / 8"                            |
| Fig-28 | 01-Pcs | DeBakey Cardiovascular Forceps 25cm / 10"                           |
| Fig-29 | 01-Pcs | Satinsky Tagential Occlusion Clamps 25cm / 10"                      |
| Fig-30 | 04-Pcs | Allis Tissues Forceps, 5x6 Teeth 25cm / 10"                         |
| Fig-31 | 01-Pcs | Russ-Modell Tissue Forceps 20cm / 8"                                |
| Fig-32 | 01-Pcs | Mayo Operating Scissors, Straight, TC Blades 17cm / 6 3/4"          |
| Fig-33 | 01-Pcs | Mayo Operating Scissors, Curved, TC Blades 17cm / 6 3/4"            |
| Fig-34 | 01-Pcs | Metzenbaum Scissors, Curved, TC Blades 20cm / 8"                    |
| Fig-35 | 02-Pcs | Metzenbaum Scissors, Curved, Delicate, TC Blades 18cm / 7"          |
| Fig-36 | 01-Pcs | Mayo-Hegar, Needle Holder, TC Jaws 18cm / 7"                        |
| Fig-37 | 01-Pcs | DeBakey, Needle Holder, TC Jaws 18cm / 7"                           |
| Fig-38 | 01-Pcs | Ryder, Needle Holder, TC Jaws                                       |
| Fig-39 | 01-Pcs | Stainless Steel Instrument Box                                      |