

Myringotomy Set

44 PCS (RAS-142-42)



Fig. 01
RC-186-07

Fig. 02
RC-186-09

Fig. 03
RD-122-03

Fig. 04
RF-126-12

Fig. 05
RF-248-12



Fig. 06
RG-586-09

Fig. 07
RG-136-12

Fig. 08
RG-160-14



Fig. 09
RI-190-16

Fig. 10
RI-262-01



Fig. 11
RI-262-03

Fig. 12
RJ-108-10

Fig. 13
RJ-126-10

Fig. 14
RJ-144-13

Fig. 15
RE-159-11

Myringotomy Set

44 PCS (RAS-142-42)



Fig. 16
RE-310-12

Fig. 17
RV-106-01

Fig. 18
RV-106-02

Fig. 19
RV-106-04

Fig. 20
RV-132-01



Fig. 21
RV-234-16



Fig. 22
RV-246-02



Fig. 23
RV-246-11



Fig. 24
RW-310-18



Fig. 25
RI-184-01



Fig. 26
RAA-104-15



Fig. 27
RF-258-10



Fig. 28
RE-123-14



Fig. 29
RE-199-14



Fig. 30
RE-206-18

Myringotomy Set

44 PCS (RAS-142-42)



Fig. 31
RL-114-13



Fig. 32
RAS-142-42
Stainless Steel Box

Myringotomy Set

44 PCS (RAS-142-42)

Sr #	Qty	Description
Fig-01	01-Pcs	Frazier Suction Tube, 2.3 mm Fr.7 19cm / 7 1/2"
Fig-02	01-Pcs	Frazier Suction Tube, 3 mm / Fr.9, Angled 19cm / 7 1/2"
Fig-03	01-Pcs	Standard Scalpal Handle, Solid, Straight No.3
Fig-04	01-Pcs	Adson Dressing Forceps, Serrated 12cm / 4 3/4"
Fig-05	01-Pcs	Adson Tissue Forceps, 1 x 2 Teeth 12cm / 4 3/4"
Fig-06	04-Pcs	Backhaus Towel Clamp Forceps 9cm / 3 1/2"
Fig-07	06-Pcs	Halsted Mosquito Forceps, Straight 12.5cm / 5"
Fig-08	02-Pcs	Kelly Forceps, Straight 14cm / 5 1/2"
Fig-09	02-Pcs	Senn Miller Retractor, Sharp, Double ended 16cm / 6 1/4"
Fig-10	01-Pcs	Langenbeck Retractor, Blade, 11x30 mm 22cm / 8 1/2"
Fig-11	01-Pcs	Langenbeck Retractor, Blade, 16x30 mm 22cm / 8 1/2"
Fig-12	01-Pcs	ALM Retractor, Sharp 10cm / 4"
Fig-13	01-Pcs	Jansen Retractor, Blunt 10cm / 4"
Fig-14	01-Pcs	Wullstein Retractor, Sharp, 3 x 3 Teeth 14cm / 5 1/2"
Fig-15	01-Pcs	Strabismus Scissors, Curved 11.5cm / 4 1/2"
Fig-16	01-Pcs	Wire Cutting Scissors, Angular, Serrated 12cm / 4 3/4"
Fig-17	01-Pcs	Hartman Ear Speculum 3 mm dia
Fig-18	01-Pcs	Hartman Ear Speculum 5 mm dia
Fig-19	01-Pcs	Hartman Ear Speculum 7 mm dia
Fig-20	01-Pcs	Ear Speculum 3 x 4 mm
Fig-21	01-Pcs	Graenwald (Jansen) Ear Forceps 16cm / 6 1/4"
Fig-22	01-Pcs	Buck Ear Curette, Fig. 00, Straight, Sharp 16.5cm / 6 1/2"
Fig-23	01-Pcs	Buck Ear Curette, Fig. 1, Straight, Sharp 16.5cm / 6 1/2"
Fig-24	01-Pcs	Freer Elevator, Double End, Sharp / Blunt 19cm / 7 1/2"
Fig-25	02-Pcs	Joseph Skin Hook, Single Prong, Sharp 15cm / 6"
Fig-26	02-Pcs	Allis Tissues Forceps, 4x5 Teeth 15cm / 6"
Fig-27	01-Pcs	Graefe Tissue Forceps, 4.5mm jaws 11cm / 4 1/2"
Fig-28	01-Pcs	Mayo Dissecting Scissors, Curved, TC Blades 14.5cm / 5 3/4"
Fig-29	01-Pcs	Metzenbaum Scissors, Curved Delicate, TC Blades 14.5cm / 5 3/4"
Fig-30	01-Pcs	Metzenbaum Scissors, Curved Delicate, TC Blades 18cm / 7"
Fig-31	01-Pcs	Halsey Needle Holder, Serrated, TC Jaws 13cm / 5"
Fig-32	01-Pcs	Stainless Steel Instrument Box