

Major Nasal Set

53 PCS (RAS-146-46)



Fig. 01
RK-162-18



Fig. 02
RF-126-12



Fig. 03
RE-218-16



Fig. 04
RE-169-15



Fig. 05
RE-150-10



Fig. 06
RE-158-11



Fig. 07
RE-159-11



Fig. 08
RW-264-16



Fig. 09
RU-116-02



Fig. 10
RW-268-15



Fig. 11
RW-310-18



Fig. 12
RW-126-35



Fig. 13
RW-284-03



Fig. 14
RW-134-50



Fig. 15
RW-284-04

Major Nasal Set

53 PCS (RAS-146-46)



Fig. 16
RW-378-01



Fig. 17
RI-184-01



Fig. 18
RW-428-12



Fig. 19
RI-184-03



Fig. 20
RW-454-16



Fig. 21
RW-450-14



Fig. 22
RW-394-10



Fig. 23
RW-460-14



Fig. 24
RW-462-16



Fig. 25
RW-482-23



Fig. 26
RW-446-11



Fig. 27
RW-348-16



Fig. 28
RV-520-04



Fig. 29
RV-234-16



Fig. 30
RW-186-16

Major Nasal Set

53 PCS (RAS-146-46)



Fig. 31
RZ-180-19



Fig. 32
RI-166-12



Fig. 33
RO-220-45



Fig. 34
RW-406-18



Fig. 35
RW-414-21



Fig. 36
RW-406-18



Fig. 37
RW-406-16



Fig. 38
RW-388-04



Fig. 39
RW-388-07



Fig. 40
RW-388-09



Fig. 41
RW-388-12



Fig. 42
RW-388-04



Fig. 43
RW-350-19



Fig. 44
RW-262-16



Fig. 45
RE-193-14



Fig. 46
RE-203-18



Fig. 47
RE-216-15



Fig. 48
RL-130-18



Fig. 49
RA-147-08



Fig. 50
RAS-146-46
Stainless Steel Box

Major Nasal Set

53 PCS (RAS-146-46)

Sr #	Qty	Description
Fig-01	01-Pcs	Farrel Cotton Applicator, Cross Serrated 18cm / 7"
Fig-02	01-Pcs	Adson Dressing Forceps, Serrated 12cm / 4 3/4"
Fig-03	01-Pcs	Cottle Nasal Scissors 16cm / 6 1/2"
Fig-04	01-Pcs	Reynolds Operating Scissors, Curved 15cm / 6"
Fig-05	01-Pcs	Iris Scissors, Straight, Sharp 10.5cm / 4 1/4"
Fig-06	01-Pcs	Strabismus Scissors, Straight 11.5cm / 4 1/2"
Fig-07	01-Pcs	Strabismus Scissors, Curved 11.5cm / 4 1/2"
Fig-08	01-Pcs	Converse Nasal Knife, Curved 16cm / 6 1/2"
Fig-09	01-Pcs	Gillies Skin Hook, Sharp 18cm / 7"
Fig-10	01-Pcs	Cottle Nasal Knife 14cm / 5 1/2"
Fig-11	01-Pcs	Freer Elevator, Double End, Sharp / Blunt 15cm / 6"
Fig-12	01-Pcs	Kilian Septum Speculum, Blades 40 mm 18cm / 7"
Fig-13	01-Pcs	Balinger Swivel Knife, Straight, 3 mm 21cm / 8 1/4"
Fig-14	02-Pcs	Cottle Septum Speculum, 50 mm Delicate Blades 15cm / 6"
Fig-15	01-Pcs	Balinger Swivel Knife, Straight, 4 mm 20cm / 8"
Fig-16	01-Pcs	Joseph Nasal Saw, Straight 18cm / 7 1/2"
Fig-17	01-Pcs	Joseph Skin Hook, Single Prong, Sharp 15cm / 6"
Fig-18	01-Pcs	Converse Alar Hooklet, Double ended, Sharp 12cm / 4 3/4"
Fig-19	02-Pcs	Joseph Skin Hook, Double prong 5 mm, Sharp 16cm / 6 1/4"
Fig-20	01-Pcs	Fomon Retractor, Two Ball End Prongs 16cm / 6 1/4"
Fig-21	01-Pcs	Cottle Double Hook, Sharp 15cm / 6"
Fig-22	02-Pcs	Cinelli Osteotome 16cm / 6 1/4"
Fig-23	01-Pcs	Cottle Alar Hooklets, 4 Prongs 14cm / 5 1/2"
Fig-24	01-Pcs	Neivert Knife Guide and Retractor, Double Ended 17cm / 6 3/4"
Fig-25	01-Pcs	Asch Septum Forceps 23cm / 9"
Fig-26	01-Pcs	Cottle Columella Forceps, 12mm wide 11cm / 4 1/4"
Fig-27	01-Pcs	Joseph Periosteal Elevator, 4 mm 16cm / 6 1/4"
Fig-28	01-Pcs	Neivert Anderson Osteotome, 4mm 20.5cm / 8"
Fig-29	01-Pcs	Gruenwald (Jansen) Ear Forceps 16cm / 6 1/4"
Fig-30	01-Pcs	Hartman Nasal Dressing Forceps 16cm / 6 1/4"
Fig-31	01-Pcs	Schmidt Haemostatic Forceps, Slight Curved 19cm / 7 1/2"
Fig-32	01-Pcs	Desmarres Lid Retractor, 12mm 13cm / 5"
Fig-33	01-Pcs	Mallet, 450 gm / 16 oz 26cm / 10 1/2"
Fig-34	01-Pcs	Maltz Rasp improved, Double ended, 8 mm wide 21cm / 8 1/4"
Fig-35	01-Pcs	Aufrecht Glabella Rasp, Single ended, 9 mm wide 21cm / 8 1/4"
Fig-36	01-Pcs	Maltz Rasp, Single ended, 10 mm wide 18cm / 7"
Fig-37	01-Pcs	Joseph Rasp, Single ended, 8mm wide 16cm / 6 1/4"
Fig-38	01-Pcs	Cottle Osteotome, 4mm 18cm / 7"
Fig-39	01-Pcs	Cottle Osteotome, 8mm 18cm / 7"
Fig-40	01-Pcs	Cottle Osteotome, 10mm 18cm / 7"
Fig-41	01-Pcs	Cottle Osteotome, 12mm 18cm / 7"
Fig-42	01-Pcs	Cottle Septum Chisel, 4mm 18cm / 7"
Fig-43	01-Pcs	Cottle Periosteal Elevator, 6mm 18cm / 7"
Fig-44	01-Pcs	Fomon Knife, Double Edge, Curved 16.5cm / 6 1/2"
Fig-45	01-Pcs	Metzenbaum Scissors, Curved, Blunt / Blunt TC Blades 14.5cm / 5 3/4"
Fig-46	01-Pcs	Metzenbaum Scissors, Curved, Blunt / Blunt TC Blades 18cm / 7"
Fig-47	01-Pcs	Fomon Dorsal Scissors, Angular, TC Blades 15cm / 6"
Fig-48	01-Pcs	Crile Wood Needle Holder, TC Jaws 18cm / 7"
Fig-49	01-Pcs	Castroviejo Caliper, 20 mm 8cm / 3 1/4"
Fig-50	01-Pcs	Stainless Steel Instrument Box