

E.R Suture Set

15 PCS (RAS-148-48)



Fig. 01
RD-122-03



Fig. 02
RF-126-12



Fig. 03
RF-248-12



Fig. 04
RC-136-12



Fig. 05
RC-138-12



Fig. 06
RC-160-14



Fig. 07
RC-162-14



Fig. 08
RC-586-09



Fig. 09
RE-104-14



Fig. 10
RE-164-12



Fig. 11
RL-130-15



Fig. 12
RAS-148-48
Stainless Steel Box

E.R Suture Set

15 PCS (RAS-148-48)

Sr #	Qty	Description
Fig-01	01-Pcs	Standard Scalpel Handle, Solid, Straight No.3
Fig-02	01-Pcs	Adson Dressing Forceps, Serrated 12cm / 4 3/4"
Fig-03	01-Pcs	Adson Tissue Forceps, 1 x 2 Teeth 12cm / 4 3/4"
Fig-04	02-Pcs	Halsted Mosquito Forceps, Straight 12.5cm / 5"
Fig-05	02-Pcs	Halsted Mosquito Forceps, Curved 12.5cm / 5"
Fig-06	01-Pcs	Kelly Forceps, Straight 14cm / 5 1/2"
Fig-07	01-Pcs	Kelly Forceps, Curved 14cm / 5 1/2"
Fig-08	02-Pcs	Backhaus Towel Clamp Forceps 9cm / 3 1/2"
Fig-09	01-Pcs	Operating Scissors, Straight, Sharp/Blunt 14.5cm / 5 1/2"
Fig-10	01-Pcs	Kilner Scissors 13cm / 5"
Fig-11	01-Pcs	Crile-Wood, Needle Holder, TC Jaws 15cm / 6"
Fig-12	01-Pcs	Stainless Steel Instrument Box