

**Overall Size (L\*W\*H):** 890\*510\*770-940mm  
**Maximum Weight Scale:** 20kg  
**Trendelenburg&Reverse Trendelenburg:** 0-10°  
**Height Adjustment:** 770-940mm

- Made of ABS Plastic Frame-Friendly Material
- Aluminum Alloy Column with Plastic Cover
- High Density and Transparent Chamber,
- Odorless, Can be Sterilized at Temperatures Above 100 °C
- Height and Tilt Function Adjustable with Gas Spring
- Rubber-Free Wheel, Wear-Resistant, Anti-Slip, Shock Absorption,
- Can Rotate 360° Freely
- Two Wheels with Brakes
- Color Option; **Pink**, **Blue**, **Green**

