

Measure the volume of bladder



KONTED



C10P

**Wireless
Bladder Scanner**

Download our ultrasound software from App Store with your ipad



3.Button for Power ON/OFF & LIVE/FREEZE

2.Battery Capacity

1.Wireless Connection Indicator

Characteristics:

- C10P Wireless Bladder Scanner is the new generation instruments for bladder volume calculating with the outstanding feature of wireless.
- The probe acts as a Wi-Fi Access Point and can be connected by iPad. With the probe be connected through Wi-Fi and the App's running, enjoy your days of working without the trouble making by cables.
- Built-in Wi-Fi and 4200mAh lithium battery



Items List for C10P Scanner

Standard Accessories:

- C10P Wireless Ultrasound Probe*1
- Portable Wireless Charger*1
- USB Cable for Charging*1
- Operators' Manual*1

Optional:

- Bluetooth printer
- ipad bracket
- ipad Trolley

1) Intended Use

Measure the volume of bladder

2) Physical characteristics

Size 180mm x 60mm x 60mm

Weight 450g

3) Electronic

Battery Capacity: 4200mAh

continuous working time: 4hour

Waterproof: IPX5

4) Probe

Frequency: 3.5MHz

Display accuracy: 1ml

Measuring range: 20 – 999ml

Measuring accuracy: 20-99ml: $\leq \pm 10\text{ml}$

100-999ml: $\leq 10\%$

5) Printer

Communication mode: Bluetooth

6) Display

Display main unit: iPad Series

Display Mode: B-Mode

Gray Scale: 256 levels