

## TECHNICAL SPECIFICATIONS

# SeaMed® LED 500

Lamp Head Diameter: 500mm

Maximum (Lighting) Light Intensity: 40.000-140.000 lux

Light Color Index: RA = 95

Color Temperature: 4300 ± 500 kelvin

Operating Power: 220/50

Lamp Energy: 60W

Lighting Setting: No Skip Setting

Focal Diameter: 160-320mm

Lamp Life (hours): 40000-50000

Lighting Depth: 700mm

# SeaMed®



## SEAMED LED-500 Surgery Lamp

DENİZLER MEDİKAL SAN. ve T. C. A. Ş.

K.bakkalköy Mh. Sevda Sk.  
Seven Tower No:1 Kat:3  
Ata ehir / STANBUL  
T : +90 216 577 52 22  
F : +90 216 577 52 21

info@denizlermedikal.com  
www.denizlermedikal.com

CE 0434

  
DENİZLER



**SeaMed**<sup>®</sup>

**LED-500**  
Surgery Lamp

### PRODUCT FEATURES

- » Lamp Power: 60W + 100W
- » Light Setting: No Skip Setting
- » Focal Diameter: 160-320mm
- » Lighting Depth: 700mm
- » Include LED500 Heads



DENİZLER MEDİKAL SAN. ve T.C.A.

K. bakkalköy Mh. Sevda Sk.  
Seven Tower No:1 Kat:3  
Ataşehir / STANBUL  
T : +90 216 577 52 22  
F : +90 216 577 52 21

info@denizlermedikal.com  
www.denizlermedikal.com

CE 0434