

TECHNICAL SPECIFICATION

SeaMed®

LED-7050

Lamp Head Diameter: 500mm + 700mm

Maximum Lighting Intensity : 140.000lux ve 170.000lux

Lighting Selection Range: %0 - %100

Light Color Index : RA =95 , RA= 95

Color Temperature : 4300±500 Kelvin

Operating Power: 220/50, 220/50

Lamp Energy: 100W, 60W

Lighting Setting: Stepless, Stepless

Focal Diameter : 160-350mm, 160-320mm

Lamp Lifetime (hour): 40000-50000, 40000-50000

Lighting Depth: 800mm, 700mm

Heads : Seamed LED700 + LED500

It prevents eye strain.

Optional camera version

DENİZLER MEDİKAL SAN. ve T. C. A. Ş.

K.bakkalköy Mh. Sevda Sk.
Seven Tower No:1 Kat:3
Ata ehir / STANBUL
T : +90 216 577 52 22
F : +90 216 577 52 21

info@denizlermedikal.com
www.denizlermedikal.com


DENİZLER

CE 0434

SeaMed®



SEAMED LED-7050

Operation Lamp



SeaMed[®]

SEAMED LED-7050
OPERATION LAMP

PRODUCT FEATURES

- » Lamp Power: 60W + 100W
- » Focal Diameter: 160-320mm + 160-350mm
- » Lighting Depth: 700mm + 800mm
- » It prevents eye strain.
- »



DENİZLER MEDİKAL SAN. ve T.C.A.Ş.

K.bakkalköy Mh.Sevda Sk.
Seven Tower No:1 Kat:3
Ataşehir / STANBUL
T : +90 216 577 52 22
F : +90 216 577 52 21

info@denizlermedikal.com
www.denizlermedikal.com

CE 0434