



GE MARQUETTE COMPATIBLE 10 LEAD BANANA ECG CABLE

Compatibility: Marquette MAC 1000, MAC 1100, MAC 1200, MAC 1600, MAC 400, MAC 500, MAC 600, MAC 800, Microsmart M, MAC 2000, MAC 3500.

Distal Connector: Screw-type DB-15 connector

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU jacket

SK1206B