



## KENZ 16-PIN COMPATIBLE 10-LEAD BANANA ECG CABLE

Compatibility: 103, 106, 107, 1201, 302, 303, Cardioline Remco Delta 1, Delta 3, Delta 30, Delta 60, Delta 1 Plus, Delta 3 Plus, Delta 60 Plus, ECG 106.

Distal Connector: Round 16-Pin Connector

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU Jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU Jacket

SK1207B