



KENZ COMPATIBLE 10 LEAD ECG CABLE

Compatibility: Kenz 109, 110, Cardico 302, 601, 1210, 1211, 108

Distal Connector: DB-15 Connector, 90° Angled

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU Jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU Jacket

SK1208B