



NIHON KOHDEN UYUMLU 10 LEAD GRABBER EKG KABLOSU

Uyumluluk: Nihon Kohden Cardiofax C, Cardiofax M, Cardiofax Q, Cardiofax S, ECG-1150, ECG-1250, ECG-1250A, ECG-1250K, ECG-1350, ECG-1350A, ECG-1350K, ECG-2150, ECG-2250, ECG-2350, ECG-9010, ECG-9020, ECG-9022, ECG-9110, ECG-9110K, ECG-9130, ECG-9130K, ECG-9132, ECG-9620, ECG-9620L.

Distal Konnektör: Vidali DB-15 konnektör

Lead Kablo Rengi: Gri

Lead Kablo Çapı: 4 mm

Lead Kablo Uzunluğu (Limb): 1.2 m

Lead Kablo Uzunluğu (V): 80 cm

Lead Kablo Malzemesi: TPU kılıf

Lead Renk Kodlaması: IEC

Lead Sayısı: 10

Proksimal Konnektör: Grabber

Trunk Kablo Uzunluğu: 3 m

Trunk Kablo Malzemesi: TPU kılıf

SK1210P ■

