



PHILIPS COMPATIBLE 10-LEAD BANANA ECG CABLE

Compatibility: Philips 200, 200i, Agilent 300i, Agilent 300PI, M1771A, M1772A, Pagewriter 100

Distal Connector: Screw-type DB-15 connector

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU jacket

SK1213B