

# SeaMed®



## SCHILLER COMPATIBLE 10-LEAD GRABBER ECG CABLE

Compatibility: Schiller AT10, AT10 Plus, AT101, AT102, AT102-Plus, AT104, AT110, AT3, AT5, AT6, AT60, C6, CS100, CS200, CS6, MS2010, MS2015

Distal Connector: Screw-type DB-15 connector

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU jacket

SK1214P