



## ZONCARE COMPATIBLE 10 LEAD BANANA ECG CABLE

Compatibility: MAC 1800, iMAC 120, iMAC 120 Pro, iMAC 12, iMAC 10, iMAC 300, E10, E12, E30, U90, ZQ-1203G, ECG-20NT

Distal Connector: DB-15 Connector with Screws

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU Jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU Jacket

SK1217B