

# SeaMed®



## MINDRAY COMPATIBLE 10 LEAD ECG CABLE

Compatibility: BeneHeart R3, BeneHeart R12, BeneHeart R12A, BeneHeart R300, BeneHeart R700, BeneHeart R900, DECG-03A

Distal Connector: Screw-type DB-15 connector

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU jacket

SK1218B