



EDAN COMPATIBLE 10 LEAD BANANA ECG CABLE

Compatibility: Edan SE-3, SE-12, SE-12 Express, SE-601A, SE-601B, SE-601C, SE-1010PC, SE-1200, SE-12 Express, ST-1212

Distal Connector: Screw-type DB-15 connector

Lead Cable Color: Gray

Lead Cable Diameter: 4 mm

Lead Cable Length (Limb): 1.2 m

Lead Cable Length (V): 80 cm

Lead Cable Material: TPU jacket

Lead Color Coding: IEC

Number of Leads: 10

Proximal Connector: Banana

Trunk Cable Length: 3 m

Trunk Cable Material: TPU jacket

SK1221B