

XM550/750

MasimoSET


Patient Monitor



Hwatime 华腾

Hwatime Medical. Make life more brilliant

The model of basic health care



XM550/XM750
Patient Monitor

Hwatime
Hwatime Medical. Make life more brilliant

7/F, Block A, Bofook Jewelry Industry Park, No.33,
Bulan Rd, Longgang District, Shenzhen City
Tel: 0086-755-28260108 Fax: 0086-755-85221277
E-mail: sales@hwatime.com



Handset two-dimensional code

Please log in for more details www.hwatime.com

Customer service hotline : 400-8972-888

· Hwatime. All Rights Reserved. The product specifications are subject to change without further notice.
· See instructions for contraindication or precautions

XM550 MasimoSET[®] Patient Monitor

Smart, Light, Easy operating
Suitable for bedside, ICU, OR



Product Feature

- Standard parameters: ECG, NIBP, RESP, PR, SpO2, TEMP
- Colorful and Clear 12.1" & 10.4" TFT screen, Backlight buttons
- Multiple display modes optional: Standard interface, Big font, ECG standard full display, OXY, Trend table, BP trend, Viewbed
- Ambulatory blood pressure technology, anti-movement
- Special design against high frequency surgical unit, and defibrillation protection
- Support Masimo SpO2-authorization partner
- 13 kinds of arrhythmia analysis
- 15 kinds of drug dose calculation
- Various languages operating systems
- Built-in detachable rechargeable lithium battery 4hours battery life
- To connect CMS, other bed observation and software updating with wireless and wired mode

Data Storage – stable and fast

- 8000 groups NIBP measurements
- 680 hours Trend data and Trend graphs
- 200 groups Alarm events reviewing
- 2 hours Wave forms reviewing
- Supporting for Arrhythmia analysis and drug dose calculation

Alarming – more safe and reliable

- 3 level Audible and visual alarming

MasimoSET[®] XM750 Patient Monitor




Movable handle design

Lightweight portable,
practical all-around,
convenient operation

XM750/550 Suitable for Adult, Pediatric and Neonatal.
Reliable quality gives confidence, stable system makes perfect.
Powerful storage capacity ensures data safety.

Hwatime Medical regards quality as the core, relies on creative R&D, produces human friendly products, and doctors' reliable partner.



NIBP list interface



Standard interface



Viewbed interface



Big font interface



Tendency diagram



OxyCRG interface



Support

HD, VGA, USB,
Internet connector output



Double invasive blood pressure (IBP)

Adopting double-channel measurement, and realizing continuous testing blood pressure variation of severe patients



• KNOB



• ONE-BUTTON NIBP



• POWER

backlight button design

Optional: EtCO2, Dual-IBP, 12-leads ECG, Touch screen, Built-in thermal printer