



DENİZLER
MEDİKAL

CATALOGUE

SeaMed LED-7050 Surgery Lamp with Camera



- * Lighting Intensity : 170.000 Lux
- * Lighting Depth: 700 mm
- * Prevents Eye Strain
- * Supported by Samsung High Quality Camera
- * 50.000 Hours Durable and Powerful LED Lamps

SeaMed LED-7050 Surgery Lamp/ Surgical Lighting



- * Lighting Intensity : 140.000 and 170.000 Lux
- * Lighting Depth: 700 mm
- * Prevent Eye Strain
- * Optional Camera Version
- * 50.000 Hours Durable And Powerful LED Lamps

SeaMed LED-500 Surgery Lamp



- * Lamp Head Diameter: 500 mm
- * Lighting Intensity : 40.000 Lux and 140.000 Lux
- * Color Index: 95 Ra
- * Electrical Power of Lamp: 60 W
- * Focal Diameter: 160-320 mm
- * Lighting Depth: 700 mm

SeaMed LED-500 Surgery Lamp



- * Lighting Intensity : 40.000 Lux and 140.000 Lux
- * Color Index: 95 Ra
- * Color Temperature: 4300 ± 500 Kelvin
- * Electrical Power of Lamp: 60 W
- * Lighting Depth: 700 mm
- * 4 Hours Battery Life

SeaMed KS-Q10 Operation Lamp



- * 78.000 Lux High Intensity Lighting
- * 10 Independent Leds
- * Color Temperature Range: 4800-5500 K
- * Easily Accessible Control Panel
- * Wall-Mounted, Mobile, Table or Ceiling Fixing Options

SeaMed HT9 Patient Monitors



- * 17" High Resolution Display (Optional Touch Screen)
- * Displaying 8-Channel Waveforms
- * Full Time Monitoring of 6 Parameters Thanks to Its Wide Screen
- * 3-Stage Audible and Visual Alarm Security
- * Wired or Wireless Connection to HIS
- * Compatible with ECG, RESP, Spo2, PR, NIBP and Temperature Measurement Probes
- * 200 Min. Internal Battery
- * Portable Design

SeaMed H9 Patient Monitors



- * 15" High Resolution Display (Optional Touch Screen)
- * Displaying 8-Channel Waveforms
- * 3-Stage Audible and Visual Alarm Security
- * Wireless Module (Optional)
- * Compatible with ECG, RESP, Spo2, PR, NIBP and Temperature Measurement Probes
- * Portable Design
- * Storing Patient Data with the aid of Internal Memory

SeaMed H8 Patient Monitors



- * 12.1" High Resolution Display (Optional Touch Screen)
- * Displaying 8-Channel Waveforms
- * 3-Stage Audible and Visual Alarm Security
- * Operating with Internal Battery
- * Wireless Module (Optional)
- * Compatible with ECG, RESP, Spo2, PR, NIBP and Temperature Measurement Probes
- * Storing Patient Data with the aid of Internal Memory

SeaMed XM550 Infant Patient Monitors



- * ECG, NIBP, RESP, PR, SPO2 and Temperature Parameters
- * Color, Sharp 10.1" or 12.1" Touch Screen Size
- * 3-Stage Audible and Visual Alarm
- * 8000 NIBP, 680 Trend Graph, 200 Group Alarm Data

SeaMed NC5 Vital Monitor



- * 8" LED Lighting LCD Screen
- * Optional Wall Mounting Solutions
- * More Than 8 Hours Working with Li-Ion Battery
- * Infant, Pediatric, Adult Mode of Use
- * Wheeled and Basket Table
- * Connection to HIS With HL7 Infrastructure



OPTOMIC OP-C5 Binocular Colposcope

- * Using with Andinclined Binocular Stereoscopic Head and Flat Binoculars
- * Halogen, Xenon and LED Light Source
- * High Lighting Intensity with Fiber Optic Cable
- * Rubber Eye Protection in Both Eyes
- * Magnification Range: 2.5x To 27x
- * Perfectly Moving Head with Angular Bearings and Teflon Discs
- * Endoscopy for Specular or Photophores and Removable Cold Light Source for Other Applications like Lightning



SeaMed Star 5000 NST Device

- * 12.1" Colored Touch Screen
- * Internal Memory
- * TOCO, FHR1, FHR FM, Internal Thermal Printer Configurations
- * 14.8 V Rechargeable Internal Battery
- * The Advantage of Connection to HIS
- * Monitoring Twin Pregnancy
- * Ability to Transfer Data to USB
- * Waterproof TOCO and Ultrasound Probes



OPTOMIC OP-C2L Binocular Colposcope

- * Using with Andinclined Binocular Stereoscopic Head and Flat Binoculars
- * Halogen, Xenon and LED Light Source
- * High Lighting Intensity with Fiber Optic Cable
- * Perfectly Moving Head with Angular Bearings and Teflon Discs
- * Endoscopy for Specular or Photophores and Removable Cold Light Source for Other Applications like Lightning



SeaMed Star 5000C NST Device

- * 12.1" Colored, Foldable LCD Touch Display
- * Monitoring ECG, Spo2, PR, NIBP and Respiratory Parameters
- * TOCO, FHR1, FHR2, FM, Internal Thermal Printer Configurations
- * 14.8 V Rechargeable Internal Battery
- * The Advantage of Connection to HIS
- * Monitoring Twin Pregnancy
- * Ability to Transfer Data to USB
- * Waterproof TOCO and Ultrasound Probes



OPTOMIC OP-C2 Binocular Colposcope

- * Using with Andinclined Binocular Stereoscopic Head and Flat Binoculars
- * Halogen, Xenon and LED Light Source
- * High Lighting Intensity with Fiber Optic Cable
- * Perfectly Moving Head with Angular Bearings and Teflon Discs
- * Endoscopy for Specular or Photophores and Removable Cold Light Source for Other Applications like Lightning



SeaMed T12 NST Device

- * 12.1" Colored LCD Touch Screen
- * Easy to Use with Functional Physical Keys
- * 3.5 Hours of Battery Life
- * Wired or Wireless Connection to HIS
- * HR, AFM, TOCO Parameters
- * Internal Probs Sockets



SeaMed KN-2200I (HD) Colposcope Device

- * Internal PC System, Printer, Trolley
- * Interface and Software are Configured for 175 Types of Cases
- * Video Recording, Image Freezing, Editing, Sectioning, Printing Reports and Storage
- * Customized Software Package Option and Compliance with Dynamic Operations
- * Zooming in, Zooming Out, Marking, Regional Growing of the Image
- * Report Editing, Printing, Prepared Templates, Diagnostic Reporting, RGB Image Printing
- * Easy Storing and Analysis with Case Number, Name, Date Informations



SeaMed KS-Q7 Examination Lamp

- * 7 LED Lighting
- * Easy Positioning with Flexible Body
- * Digital Brightness Control
- * Head That Can Move 360°
- * Whelled Mobile Stand

SeaMed KN-2200A Colposcope

- * Compatibility with Win7 and Win8 Operating System
- * Video Recording, Image Freezing
- * Regulation, Sectioning, Report Printing and Storing
- * Zooming in, Zooming Out, Marking, Regional Growing of the Image
- * Report Editing, Print, Ready Templates
- * Easy Storing and Analysis with Case Number, Name, Date Informations



SeaMed KS-Q5 Examination Lamp

- * 40 mm Lamp Head Diameter
- * Gradual Light Adjustment With Dimmer
- * 1 LED Light Source
- * Angle Adjustment On The Lamp Head
- * Whelled Mobile Stand



SeaMed Headlamp

- * Internal Rechargeable Lithium Polymer Battery
- * Adjustable Spot Size
- * Optional Yellow Filter
- * Digital Brightness Control

SeaMed Fetal Doppler



- * Compact And Portable Design
- * Wide Lighting LCD Display
- * Replaceable Waterproof Prop
- * Internal Speaker
- * Full Charged Battery Operating More Than 10 Hours Continuously
- * High Sensitive

SeaMed Water Birth Pool



- * Large and Medium Size
- * 6 Handles
- * High Strength and Healty Material in 6P-PVC Quality
- * Including Water and Air Pump
- * Suitable for Professional Use

Seamed BS6930 ECG Holter



- * High Accuracy
- * 7 Days Recording Capacity
- * Comfortable and Lightweight Design
- * Suitable for Operating with both 12 and 3 Channel
- * Sd Card Supported
- * Easy Battery Replacement with Alkaline Battery

SeaMed iTengo + ECG Holter Software



- * Multi-channel Operation
- * % 99 Consistency
- * Displaying PR Ranges with High Accuracy (uV)
- * Detecting Small Pathological Changes
- * T-Wave Period Analysis
- * Analysis Tracking in Multiple Window
- * ST Segment Analysis
- * R-R Trend and Its Histograms
- * Waterfall Analysis
- * ASCII ve GDT Data Support



Seamed Echo Blood Pressure Holter

- * Comfortable and Silent Operation
- * Long Battery Life due to Low Power Consumption
- * Posture Detection via Sensors (Lying, Sitting, Walking)
- * Carry Case and Sleeves Suitable for Using
- * Easy Analysis With Software
- * Ergonomic Structure For Daily Use

SeaMed H12 ECG Device



- * 10.4" LED HD Touch Screen
- * Automatic Arrhythmia Detection and Report
- * Reporting DICOM, BMP, JPG, PDF, SCP-ECG, FDA-XML File Format
- * Integration into the HIS
- * Entering Patient Data by Handwriting
- * Monitoring and Printing 600 Second ECG Waveforms
- * Pace-maker Detection
- * Easy to Use with Touch Pen

SeaMed H3 ECG Device



- * 6.2" LCD Touch Screen
- * 3x4 ve 6x2 Screen Format
- * Integration into the HIS
- * Automatic Arrhythmia Detection
- * 122 Type Diagnosis Report
- * Pace-maker Detection
- * Monitoring and Printing 600 Second ECG Waveforms
- * Internal WIFI
- * Operating More Than 10 Hours Continuously with Full Battery

SeaMed 1200B ECG Device



- * 5.7" TFT Screen
- * Simple and Clear Keyboard Order, Backlight Support
- * Internal Memory with 300 ECG Recording
- * Digital Filter with High Accuracy
- * External Printer Connection and USB Port for Data Transfer
- * PC-ECG Management System

SeaMed 300B ECG Device



- * 3.5" Colored LCD Screen
- * Internal Thermal Printer
- * Internal Rechargeable Li-ion Battery
- * Digital Filter
- * Internal Memory with 300 ECG Recording
- * External Printer Connection and USB Port for Data Transfer
- * 1/2/3 Channel ECG Waveform

SeaMed Capnography Device



- * Using in Infant, Child and Adult Patients
- * Monitoring EtCO₂ / CO₂, Respiration Parameters
- * Monitoring EtCO₂ Graphically
- * Measuring in 10 Seconds
- * Audio and Visual Alarm

SeaMed KL-20 Oxygen Air Blender



- * Compact Lightweight Design
- * Adjustment of Oxygen Rate with One Rotary Button
- * Operating Between 21% and 100% Rates
- * Flow Adjustment Between 0 L / Min. and 15 L / Min.
- * Operating Between 2,8 and 6 Bar
- * Fixation Vertically with Apparatus
- * Using with Infant Incubator

SeaMed PC900B Capnography Device



- * Using in Adult, Child and Infant Patients
- * Monitoring EtCO₂ / CO₂, Respiration, Spo₂, PR Parameters
- * Monitoring EtCO₂ and Spo₂ Parameters Graphically
- * Using in Intubated and Non-Intubated Patients
- * Audio and Visual Alarm
- * Data Transfer to Computer via USB Port

Seamed FT2800 Dual Blood and Infusion Warmer



- * High Performance Dry Heating System
- * Microprocessor Control
- * Standard Set for Medical Fluids
- * With Heater Profile Holder
- * Operating in High Accuracy with the aid of Integrated Temperature Controller
- * Ease of Operation with Dual-Channel
- * Easy Tracking with the aid of Wide LCD Screen

SeaMed LA204 Hand Type Spirometer



- * Easily Tracking with the aid of LCD Colored Screen
- * Rechargeable Lithium Battery
- * 0-16L/S Flow Scale
- * Increased Accuracy With Digital Turbine Sensor
- * Connecting to P hone With Bluetooth
- * 315 Gram Lightweight and Compact Design

Seamed FT1800 Blood and Infusion Warmer



- * Using in Adult, Child and Newborn Patients
- * Hızlı Isınma ve Isıtma Sistemi
- * Temperature Adjustment Between 37° C and 42° C
- * Overheating Security System
- * Possibility of Fixing Serum Hanger
- * Easy Tracking with the aid of Wide LCD Screen
- * High Accuracy Operating with Integrated Temperature Control

SeaMed Pc Based Spirometer



- * Accuracy of Test Results Comply with ATS Standards
- * Standard Estimates for Different Ethnic Groups and Lung Age Decision
- * Built-in Bronchial Expansion and Provocation Test Modules
- * Live Help Animation to Improve Test Collaboration

Seamed FT800 Blood and Infusion Warmer



- * Effective Operating with 3 Temperature Sensors and Microprocessor Controlled, High Security
- * Fast, Automatic Use with One Key
- * Compatible with Standard Infusion and Blood Sets
- * Easy, Fast Installation and Using
- * High Accuracy Operating with Integrated Temperature Control

SeaMed LA104 Console Type Spirometry



- * 8.4" Touch Screen
- * Durable Spiral Cable Connecting the Probe to Console
- * Easy Calibration
- * Faster Measurement with Dual Processor
- * Ability to Store 8000 Patient Recordings
- * USB Card Supported

SinoMDT SN-H10 Blood And Infusion Warmer



- * High Performance Dry Heating System
- * Microprocessor Controlled
- * Standardized Sets for Medical Fluids
- * Sensor-Controlled Electrical Heating Profiles
- * Safe Use with Visual and Audible Alarm
- * Interface for Monitoring Temperature, Warm-Up Time, Patient Status
- * Suspension Bracket for Serum

SeaMed Anim-8 Infusion Pump



- * 3.5" Colored TFT, LCD Touch Screen
- * Set Independent Horizontal Peristaltic Pump
- * Ability to Operate With All Infusion Sets
- * Turkish User Menu
- * Microprocessor Controlled Reliable and Precise Performance
- * KVO, Bolus and Purge Functions
- * Audio and Visual Alarm
- * Optional Docking Station System
- * Connection to HIS

SeaMed Anim-5 Injector Pump



- * 3.5" Colored TFT, LCD Touch Screen
- * Set Independent Structure
- * Turkish User Menu
- * 2.5, 10, 20, 30, 50/60 Volume Injector Usage
- * Microprocessor Controlled Reliable and Precise Performance
- * Audio and Visual Alarm
- * Optional Docking Station System
- * Connection to HIS

SeaMed Semi-2 Infusion Pump



- * 3.5" Colored TFT, LCD Touch Screen
- * Set Independent Horizontal Peristaltic Pump
- * Ability to Operate With All Infusion Sets
- * Turkish User Menu
- * Microprocessor Controlled Reliable and Precise Performance
- * KVO, Bolus and Purge Functions
- * Audio and Visual Alarm
- * Optional Docking Station System

SeaMed Semi-1 Injector Pump



- * 3.5" Colored TFT, LCD Touch Screen
- * Set Independent Form
- * Turkish User Menu
- * 2.5, 10, 20, 30, 50/60 Volume Injector Usage
- * Microprocessor Controlled Reliable and Precise Performance
- * Audio and Visual Alarm
- * Optional Docking Station System

Zerone INFU-Z 2000 Dual Infusion Pump



- * Touch Keypad
- * Ergonomic Design
- * Quick and Easy Installation
- * Suitable Weight For Easy Carrying
- * KVO and Bolus Functions
- * Microprocessor Controlled Reliable and Precise Performance
- * Audio and Visual Alarm
- * Easy Battery Replacement
- * Silent Operation

SeaMed KL-602 Injector Pump



- * 2.3" Colored TFT, LCD Display
- * Initializing, Flow Rate, Temporal, Body Weight, Modes Of Use
- * Set Independent Form
- * Turkish User Menu
- * 10, 20, 30, 50/60 Volume Injector Usage
- * Microprocessor Controlled Reliable and Precise Performance
- * KVO, Bolus and Purge Functions
- * Audio and Visual Alarm

SeaMed KL-8052 Infusion Pump



- * 4" Colored TFT, LCD Display
- * Flow Rate, Drop and Temporal Usage Modes
- * Set Independent Peristaltic Pump
- * Ability to Operate with All Infusion Sets
- * Turkish User Menu
- * Microprocessor Controlled Reliable and Precise Performance
- * Audio and Visual Alarm

SeaMed KL-702 Injector Pump



- * 3" Colored TFT, LCD Display
- * Microprocessor Controlled Reliable and Precise Performance
- * Initializing, Flow Rate, Temporal, Body Weight, Modes Of Use
- * Set Independent Form
- * 10, 20, 30, 50/60 Volume Injector Usage
- * Dual Channel Transmission at Different Speeds
- * KVO, Bolus and Purge Functions
- * Audio and Visual Alarm

SinoMDT SN-600N Feeding Pump



- * Convenient For Cleaning Thanks to Its Transparent Cover
- * Safe Use Thanks to Its Waterproof Feature
- * Adjustable Gradual Audible and Visual Alarm System
- * Easy Calibration Thanks to DUAL-CPU Feature
- * USB Support

Ampall PCA Pump



- * Microprocessor Controlled Reliable and Precise Performance
- * Turkish User Menu
- * Save Up to the Last 15 Operation
- * Recording Drug Name and Infusion Information
- * 3 Stages Air and Occlusion Level
- * KVO, Bolus and Purge Functions
- * Audio and Visual Alarm
- * 50 Hours of Battery Life



Cheiron Victoria Lipos Liposuction Device

- * Silent and non-vibrating Operating in Liposuction Operations
- * Durable and UV Protected Design
- * Easy to Use with Footswitch
- * 39.4 dB Sound Level
- * Microbiological Bacteria Filter
- * Maintenance-Free Oil-Free Pump
- * Suitable Design for Operating Continuously
- * Cooling With Extra Fan Against Heating



SeaMed 7A23-B Surgical Aspirator

- * Easy to Use with Integrated Footswitch
- * Easy Transition Between Jars
- * Maintenance-Free Oil-Free Pump
- * Suitable Design for Operating Continuously
- * Easy to Carry
- * Low Noise



Cheiron Victoria Thorax Drainage

- * Durable and UV Protected Design
- * 39.4 dB Sound Level
- * Microbiological Bacteria Filter
- * Maintenance-Free Oil-Free Pump
- * Suitable Design for Intermittent Operation
- * Cooling With Extra Fan Against Heating



SeaMed 7E-G Transport Aspirator

- * Continuously Operating with Intermittent Loading
- * Silent And Vibration-Free Operating in Mini Aspiration Processes
- * Compact Lightweight Transport
- * Transport Protection System
- * Maintenance-Free Oil-Free Pump
- * Suitable Design for Operating Continuously



Cheiron Victoria Versa Surgical Aspirator

- * Silent and Vibration-Free Operation in All Surgical Procedures
- * Durable and UV Protected Design
- * Easy to Use with Footswitch
- * 39.4 dB Sound Level
- * Microbiological Bacteria Filter
- * Maintenance-Free Oil-Free Pump
- * Suitable Design for Operating Continuously
- * Cooling With Extra Fan Against Heating



Dr. Oppel ST-501 Radiofrequency Device

- * Precise Cutting and Coagulation with 4 Mhz Radio Frequency Current
- * Using In ENT, General Surgery, Dermatology, Urology, Plastic Surgery
- * Minimal Damage in Environmental Tissues
- * Time Selection: Continuous, 1 sec, 1/2 sec, 1/4 sec, 1/8 sec
- * Monopolar Operating Power: 130 W
- * Bipolar Operating Power: 90 W

Ca-Mi Hospivac Surgical Aspirator

- * Silent and Vibration-Free Operation in All Surgical Procedures
- * Easy to Use with Footswitch
- * Microbiological Bacteria Filter
- * Maintenance-Free Oil-Free Pump
- * Suitable Design for Operating Continuously
- * Protection Against Heating Problem



Edan IT-20 Telemetry

- * Lightweight Compact Structure
- * Easy-to-Read OLED Screen
- * Monitoring ECG, Respiration, Spo2, PR Parameters
- * 3 or 5 Lead ECG Cable
- * Different Display Mode Options
- * Connecting to System with 64 Devices at the same time
- * Nurse Call and Patient Call Functions
- * Connecting to Central System with WIFI
- * Recording 240 Hours Trend and 720 Hours Alarm History for Each Patient



Ca-Mi Askır 36 Br Transport Aspirator

- * Silent And Vibration-Free Operation in Mini Aspiration Processes
- * Transport Usage with Compact Lightweight
- * 3.5 Hours Operating With Internal Battery
- * Maintenance-Free Oil-Free Pump
- * Suitable Design for Continuous Operation



Bistos BT-200 Vascular Doppler

- * Lightweight Compact Structure and Easy Operation
- * 4, 5 And 8 Mhz Probes for Different Clinical Applications
- * High Precision Probes
- * Internal Headphone Output
- * Diabetic Foot Application



UTAS Uvent A Intensive Care Ventilator



- * High Resolution Colored Touch Screen
- * Monitoring Vital Symptom like Spo2 ve CO2
- * Audio and Visual Alarm
- * Various Ventilation Modes
- * Unsaturated Oxygen Sensor
- * Modes: VCV (S), PCV (S), V-SIMV, PSIMV, BiPPV/APRV, CPAP, nCPAP, PS, PCVT/PSVT, APNEA

Seamed DH-320 Humidifer



- * Visual and Adjustable Limits Alarm
- * Compact Design with Colored Wide LCD Screen
- * Dual Temperature Control Feature
- * Unique PID Control System
- * Ability to Operate Without Accumulating Water in The Circuit
- * Ability to Reach Quickly to Target Temperature
- * Dual and Single Limb Option
- * Invasive and Noninvasive Usage Option

UTAS Uvent T Intensive Care Ventilator



- * Classic and Smart Ventilation Modes for Infant, Child and Adult Patients
- * Various Ventilation Modes
- * Internal Additional Monitoring Channels
- * Various Ventilation Modes
- * Unsaturated Oxygen Sensor
- * Compact Design
- * Usigng with VCV, CPAP, V-SIMV, PCV, PCVT, Ncpap/HFOT Modes

Seamed SN-100 Humidifer



- * Electronic Chip Controlled Temperature Tracking
- * Protection System Against Water Spills
- * Chamber Locking Mechanism
- * Heating Condition LED Indicator
- * Secondary Security System with Physical Thermostat
- * Adjustable Thermostat Adjustment Between 28-37 Degrees

UTAS Uvent M Transport Ventilator



- * Microprocessor Controlled Reliable Ventilation
- * Smart Functional 12" Colored Led Touch Screen
- * Registration up to 10,000 Patients
- * Ergonomic Using with Standby, Lock, Freze Modes
- * Sensitivity Range with Wide Breath Range

Insighters IS3-L Video Laryngoscope



- * 3.5" Led Screen
- * LED Light Source
- * ANTI-FOG Feature
- * Durable Optical Polymer Blades
- * Anatomical Structure Does Not Pressure the Ideal Teeth
- * Easy to Use And Practical
- * Disposable Blades

Seamed Oxi 4 Transport Ventilator



- * Smart Functional Touch 7.1" LCD Screen
- * Suitable for Adult and Pediatric Patients
- * Suitable for Patient Transfer with Ergonomics
- * Simple and Innovative, Safe Design
- * FIO2 Regulation Range Between 40 % -100 %

Insighters IS3-R Rigid Video Laryngoscope



- * 3.5" Led Screen
- * LED Light Source
- * Vapor Preventer ANTI-FOG Feature
- * Anatomical Structure Does Not Pressure the Ideal Teeth
- * Easy to Use And Practical

Seamed Oxi 2 Plus Transport Ventilator



- * 5" Colored Touch Screen
- * Suitable for Adult and Pediatric Patients
- * Simple, Innovative and Ergonomic Design
- * FIO2 Regulation Range Between 50% And 100%

Insighters IS3-F Flexible Video Laryngoscope



- * Using in Infant, Child and Adult Patients
- * Data Transfer with USB and Charging with USB Socket Input
- * Image Transfer to External Monitor with HDMI
- * Optional Wi-Fi
- * Record Storage with 8 GB Memory
- * Retrospective Analysis with Recording Property
- * Easy and Fast Installation



Heal Force 8502D Infant Incubator

- * Internal Weighing and Phototherapy Unit
- * Kejian KJ-8000 Bilirubinmeter
- * New Generation Fiber, Photoelectric and Data Processing Technology
- * Getting Quick Results Without Disturbing The Baby
- * Patented Aerodynamic Design
- * Easy-to-Read LCD Screen
- * Full Battery Charge in 3 Hours With Quick Charge
- * Low Battery Warning
- * Accuracy Equivalent to High Performance HPLC Method



Kejian KJ-8000 Bilirubinmeter

- * New Generation Fiber, Photoelectrical and Data Processing Technology
- * Quick Results Without Disturbing The Baby
- * Patented Aerodynamic Design
- * Easy-to-Read LCD Screen
- * Full Battery Charge in 3 Hours with Fast Charge
- * Low Battery Warning Icon
- * Accuracy Equivalent to High Performance HPLC Method
- * Automatic Calibration Feature



Heal Force BN-3000A Infant Incubator

- * Servo Control Air Flo
- * Temperature Indicator
- * Heater Power Indicator
- * Audible and Visual Alarm
- * Pre-installation Mode and Self-test Feature
- * Optional Skin Temperature and Humidity Indicator



SonoScape P20 Ultrasound

- * Super Band width Platform
- * Spatial Composite Imaging
- * μ -Scan
- * 4D Diagnostic Technology
- * Touch Screen
- * Application Field: General, Radiology, Cardiology, Obstetric and Gynecology, Urology Vascular



Heal Force FXQ3 Radiant Warmer

- * Preliminary Function With APGAR Timer
- * Swinging Bearing and Rotating Radiant Head
- * X-Ray Cassette Mount
- * Heating Feature With Anti-Explosion Quartz Infrared Tube
- * Self-Test Feature
- * Visual and Audio Alarm



SonoScape X3 Portable Ultrasound

- * Small, Portable Ergonomic Design
- * μ -Scan, Multiple Beam Processing, Spatial Composite Imaging, Pure Inversion, Harmonic Imaging
- * Using with Convex, Linear, Echo, Intracavitary and Microconvex Probes
- * Application Fields: General, Radiology, Cardiology, Obstetric and Gynecology, Urology Vascular



Heal Force GLQ01 Led Phototherapy

- * Blue Tube or LED Compatible Therapy Supply
- * LED Digital Display
- * Digital Therapy Time Indicator
- * Rotatable Radian Head 30/60/90 Degrees



SeaMed Convex HandHeld Ultrasound Device

- * Black&White Ultrasound Featured
- * B and B/M Visualization
- * 80 Ceramic Piezo Electrical Elements
- * 3.5 / 5 Mhz, R60 Convex Graded Probe
- * 260 Gram Lightweight and Compact Design



SeaMed Tunnel Type Phototherapy

- * 5 "Touch LCD Screen Control Panel
- * Equal Light Distribution to Every Point
- * 4 Levels Light Intensity Setting
- * Adjustable Alarm Control Unit
- * Continuous Ambient Temperature Measurement
- * Baby Skin Temperature Measurement
- * Easy Access to Baby



SeaMed Linear Prob Hand Type Ultrasound Device

- * Black&White Ultrasound Featured
- * B and B/M Visualization
- * 80 Ceramic Piezo Electrical Elements
- * 7.5 / 10 Mhz, L40 Convex Graded Probe
- * 260 Gram Lightweight and Compact Design



SeaMed S8 Defibrillator

- * 8.4 " Colored TFT Screen
- * Displaying Maximum 2 Waveform
- * 800 × 600 Pixels
- * Two Rechargeable Li-Ion Batteries
- * Manual Defibrillation, Synchronized Defibrillation and AED Modes
- * 200j Charging in less than 5 Seconds, 360j Charging in less than 8 Seconds



SeaMed Wrist Sphygmomanometer

- * Smart Pressure Control Technology
- * Wide LCD Screen
- * Automatic Inflating and Deflating
- * Powered by AA Size Alkaline Battery
- * Automatic Shutdown Function



SeaMed CM-1329 Wall Mounted Sphygmomanometer

- * Wall-Mountable Aneroid Type
- * Easily Readable Display
- * 0-300 mmHg Measuring Scale
- * Easy Mounting Bracket with Top Slide and Plug-in Type
- * Spiral Hose That Can Extend Up to 3 Meters
- * Healthy Design Without Mercury
- * With Air Discharge Tap



SeaMed CM-4236 Pediatric Stethoscope

- * Diaphragm Made of Stainless Steel
- * Extra Large Bell for Frequency Response
- * Ultra Sensitive Diaphragm for Clear High Frequency Response
- * Superior Wearing Comfort and Acoustic Seal
- * PVC Ear Tips



SeaMed CM-3003 Sphygmomanometer

- * Nickel-Plated and Chrome-Plated Manometer Chamber
- * 0-300 mmHg Measuring Scale
- * Dual Tube Technology
- * Luer Connector System
- * Easily Inflatable Flexible Cuff Axis
- * Safety, Fixed Air Drain Screw
- * Chrome Plated Brass Relief Valve
- * Comes with Carrying Case and Accessories



SeaMed CM-4136 Adult Stethoscope

- * Diaphragm Made of Stainless Steel
- * Extra Large Bell for Frequency Response
- * Ultra Sensitive Diaphragm for Clear High Frequency Response
- * Superior Wearing Comfort and Acoustic Seal
- * PVC Ear Tips



SeaMed CM-3003 Pediatric Sphygmomanometer

- * Nickel-Plated and Chrome-Plated Manometer Chamber
- * 0-300 mmHg Measuring Scale
- * Luer Connector System
- * Safety, Fixed Air Drain Screw
- * Chrome Plated Brass Relief Valve
- * Comes with Carrying Case and Accessories



SeaMed CM-4113 Adult Stethoscope

- * Extra Large Bell for Frequency Response
- * Ultra Sensitive Diaphragm for Clear High Frequency Response
- * Superior Wearing Comfort and Acoustic Seal
- * PVC Ear Tips



SeaMed Upper Arm Digital Sphygmomanometer

- * Smart Pressure Control Technology
- * Wide LCD Screen
- * Automatic Inflating and Deflating
- * Powered by AA Size Alkaline Battery
- * Automatic Shutdown Function



SeaMed Forehead Thermometer

- * Operating with 2x AAA Battery
- * Body And Object Measurement Modes
- * LCD Color Warning by Temperature
- * Aural Warning
- * Lightweight and Compact Design

Covidien Console Type Pulse Oximeter



- * Easy, Intuitive Operation and Small Design
- * 10-Hour Battery Life
- * 96-Hour Trend Memory Recorded Every Four Seconds
- * Internal Memory Recording
- * Self-Holding Handle, Small, Portable, Lightweight, Easy to Carry Design
- * Multi Language Support Graphical User Interface



SeaMed Fiber Optic Macintosh Laryngoscope

- * Safety Use in Patients
- * Easy Xenon LED and Battery Replacement
- * 3500 (2.5 V) or 10000 (3.5 V) Lux Light Intensity
- * Non-sharp Curved Blades
- * 4000 Times Octavable At 134° C Temperature
- * 5 Type Macintosh Blades
- * Stainless Steel Structure

SeaMed SP-20 Handheld Pulse Oximeter



- * Colored TFT LCD Display
- * Usage in Adult, Pediatric, Newborn Patients
- * Fast and Reliable Measurement
- * Audio and Visual Alarms
- * Different Screen Display Mode
- * Optional Fever Meter



SeaMed Fiber Optic Miller Laryngoscope

- * Safety Use in Patients
- * Easy Xenon LED and Battery Replacement
- * 3500 (2.5 V) or 10000 (3.5 V) Lux Light Intensity
- * Non-sharp Curved Blades
- * 4000 Times Octavable At 134° C Temperature
- * 5 Type Macintosh Blades
- * Stainless Steel Structure

SeaMed GM100P20 Finger Type Pulse Oximeter



- * Colored OLED Display
- * Use in Adult, Pediatric Patients
- * Fast and Reliable Measurement
- * Audio and Visual Alarms
- * 4 Different Screen Display Modes
- * Automatic Power On / Off
- * Up to 50 Hours of Study



SeaMed Maxilite Macintosh Laryngoscope

- * Safe Use in Patients
- * Closed Fiberoptic system
- * Easy Battery and Xenon LED Replacement
- * 3500 (2.5 V) or 10000 Lux (3.5V) Light Intensity
- * Non-Sharp Curved Blades
- * 4000 Times Octavable At 134° C Temperature
- * 5 Type Macintosh Blades
- * Stainless Steel Structure

SeaMed POD4 Pediatric Pulse Oximeter



- * Lightweight, Small, Easy to Carry with Carrying Case and Strap.
- * Accurate Measurement of SpO2 and Pulse Rate Values;
- * SpO2, Pulse Rate, PI, Heart Rate and Plethysmogram Display on the Colored Screen
- * 4-Directional Screen Display and Various Display Modes;
- * Automatic Measurement Start After Insertion
- * Automatic Signal Off If No Signal For 4 Seconds
- * More than 30 Hours of Continuous Working Time with 2 AAA Batteries
- * Visual and Audio Alarm



SeaMed Maxilite Miller Laryngoscope

- * Safe Use in Patients
- * Closed Fiberoptic system
- * Easy Battery and Xenon LED Replacement
- * 3500 (2.5 V) or 10000 Lux (3.5V) Light Intensity
- * Non-Sharp Curved Blades
- * 4000 Times Octavable At 134° C Temperature
- * 5 Type Miller Blades
- * Stainless Steel Structure

SeaMed Hard Intubation Laryngoscope Set



- * Handy Type for Difficult Intubations
- * Flexi Type Macintosh Blade That Can Bend 70 Degrees
- * Prevention of Possible Tooth and Soft Tissue Destruction
- * Integrated Fibers for Best Lighting
- * High Quality Stainless Steel Structure
- * Autoclave can be Entered



SeaMed SM-140 Otoscope Set

- * Optical Lighting with LED Technology
- * Metal Compact Durable Design
- * 3 x Magnifying Durable Window
- * Adjustable Light Intensity
- * Special Protection Bag

SeaMed SM-147 Ophtalmoscope Set



- * Lighting with LED Technology
- * Metal Compact Durable Design
- * -20 +20 Diopter Lens
- * Integrated On/Off Switch
- * Adjustable Light Intensity
- * Special Protection Bag

SeaMed Wall Mounted Otoscope Set



- * Automatic Light-on When Lifted
- * 10 Stages Adjustable Light Intensity (18.5 Lumens)
- * Ergonomic Handling Thanks to the Extended Spiral Cable
- * Durable Outer Plastic Housing
- * Easy-to-Remove Speculum
- * Internal Speculum Chamber

SeaMed Otoscope Set



- * Remote Fiber Optic Lighting with LED Technology
- * Compact Durable Design
- * 3 x Magnifying Durable Window
- * Speculum in 4-D
- * Special Protection Bag

SeaMed Wall Mounted Ophtalmoscope Set



- * Automatic Light-on When Lifted
- * 10 Stages Adjustable Light Intensity (18.5 Lumens)
- * Blue, Green and Natural Light Filter
- * -20 +20 Diopter Range
- * Ergonomic Handling Thanks to the Extended Spiral cable
- * Durable Outer Plastic Housing

SeaMed Ophtalmoscope Set



- * Lighting with LED Technology
- * Compact Durable Design
- * -20 +20 Diopter Lens
- * Integrated On/Off Switch
- * Special Protection Bag

SeaMed Air Pressure Mattress (Rotational)



- * Effective Treatment That Covers The Whole Body
- * A + B + C Operating System
- * Integrated CPR Module
- * Possibility to Use Up to 200 Kg
- * Easy Connection and Installation
- * 10/15/20/25 Min. Cycle Time
- * Waterproff Air-Permeable Antibacterial Cover

SeaMed SM-173 ENT Diagnostic Set



- * In Special Box
- * Octavable at 134° C
- * Inside Adjustable Lighting Set
- * 1 Otoscope
- * 1 Ophtalmoscope
- * 2 Laryngeal Mirror(3 and 4 Sizes)
- * 3 Ear Speculum (2,5 / 3,5 / 4,5 Mm Size)
- * 1 Light Holder (For Laryngeal Mirror)
- * 1 Nose Holder (Reusable)

SeaMed Air Pressure Mattress Tubular



- * Effective Treatment That Covers The Whole Body
- * A + B + C Operating System
- * Integrated CPR Module
- * Easy Connection and Installation
- * Possibility to Use Up to 150 Kg

SeaMed Wall Mounted Otoscope-Ophtalmoscope Set



- * Automatic Light-On When Lifted
- * 10 Stages Adjustable Light Intensity (18.5 Lumens)
- * Ergonomic Using Thanks to the Extended Spiral Cable
- * Durable Outer Plastic Housing
- * Easy-to-Remove Speculum
- * Internal Speculum Chamber
- * Blue, Green and Natural Light Filter
- * -20 +20 Diopter Range

SeaMed Air Pressure Mattress Balloon Type



- * Bubble Type PVC Bed And Pump
- * Powerful Pump Unit
- * Adjustable Pressure Operability
- * Can Be Used 24 Hours A Day
- * Provides Full Support to Patients Weighing Up to100 Kg

DENİZLER MEDİKAL SAN. ve TİC. A.Ş

**K.Bakkalköy Mah. Sevda Sk.
Seven Tower No:1 K:3
Ataşehir / İSTANBUL**

Tel: +90 216 577 52 22 - 25

Fax: +90 216 577 52 21

info@denizlermedikal.com

www.denizlermedikal.com